

Spin Selective Neutron Transportation and Shaping by using Multiplet Magnetic Lenses

JAEA : T. Oku, J. Suzuki, H. Kira, T. Shinohara, S. Takata,
M. Arai

RIKEN : K. Hirota

Kyoto Univ. : Y. Iwashita

KEK : H.M. Shimizu



Contents:

- 1) Neutron Polarizing Guide based on the Quadrupole Magnet.**
- 2) Pulsed Neutron Beam Focusing by using Multiplet-Magnetic Lens System for TOF SANS instruments.**



Neutron Motion in a magnetic field

Potential Energy : H

$$H = -\mu \cdot B \quad \rightarrow \quad m \frac{d^2 r}{dt^2} = -\nabla \mu \cdot B$$

Magnetic field torque neutron spin


$$\frac{d\sigma_n}{dt} = \gamma \sigma_n \times B, \quad \sigma_n = \mu / |\mu|$$



Larmore Precession

Larmore Angular Frequency : $\omega_L = \gamma |B|$

Angular Frequency of magnetic field vector at neutron position : $\omega_B = \frac{ds}{dt} \left| \frac{\partial}{\partial s} \frac{B}{|B|} \right|$



If the magnetic field is strong enough,

$$\frac{\omega_L}{\omega_B} \gg 1$$



Spin component parallel to the magnetic field is almost conserved


Positive polarity

$$\frac{d^2 r}{dt^2} = + \frac{\mu}{m} |\nabla |B||$$

Negative polarity

positive polarity negative polarity

B spin B spin




Neutron spin polarity




Neutron motion in magnetic field

Condition :
 $|IB|=2T$
 $@r=2.5\text{mm}$





Neutron motion in magnetic field

Condition :
 $|IB|=2T$
 $@r=2.5\text{mm}$



Polarizing power

condition : $|IB| = 2T @ r = 2.5\text{mm}$




Polarization degree P as a function of the Elapsed time.
Here, P was calculated by counting the neutron numbers with up
and down spins within $r < 4\text{mm}$ after the elapsed time.



Halbach-type Quadrupole Magnet


Cd cylinder inserted into the magnet bore




$$IBI = G_q r, \quad G = 800 \text{ T/m}$$



Acceleration of neutron
 $\sim 4600 \text{ m/s}^2 \text{ } (\sim 470G)$




Weight
 $\sim 25\text{kg}$




Cold Neutron Polarization Experiment

Experimental setup




Neutron Polarization inside the Quadrupole Magnet




*Magnetic field distribution of
quadrupole magnet*

Neutron Polarization inside the Quadrupole Magnet




Cold Neutron Polarization Experiment

Experimental setup




Polarization Analysis using the Magnetic Lens

Negative Polarity



Negative Polarity




Stern-Gerlach experiment

Polarization degree :

$$P = \frac{I_+ - I_-}{I_+ + I_-}$$



Polarization Analysis using the Magnetic Lens




Stern-Gerlach experiment

$$\text{Polarization degree : } P = \frac{I_+ - I_-}{I_+ + I_-}$$

Experimental Result


Setup

Beam line : JRR-3 C3-1-2-1 NOP




off

on




Plan of Practical Application

The Quadrupole Magnet are planed to be installed at the Focusing SANS Instrument SANS-J-II in March 2007 by JAEA Neutron Optics and Soft-matter Groups.




Halbach-type quadrupole magnet



Focusing SANS Instrument
with the Magnetic Lens
(SANS-J-II of JRR3)

Neutron Polarizing and Transporting Devices of SANS-J-II




Ni-Guide
Polarizing Mirror 1
Polarizing Mirror 2
Quadrupole Magnet

Variable-Size Slit (Slit1)

Inside of the Collimator chamber of SANS-J-II


Comparison

Magnetic mirror vs Quadrupole magnet





Characteristic of the Quadrupole Magnet Polarizer

- ❖ *High polarizing power \(^{\wedge}\circ\wedge)*
 - ❖ *High Polarization degree of $P=0.9993\pm0.0025$ was achieved.*
- ❖ *High stability o(^{-}\wedge)o*
 - ❖ *It is composed of the permanent magnets.*
- ❖ *High transmission d(\wedge o\wedge)b*
 - ❖ *No attenuation of polarized neutrons occurs, since they do not interact with any materials.*
- ❖ *Opening aperture is small. (T.T)*
 - ❖ *The magnet aperture can be enlarged to some extent.*
- ❖ *Magnet length must be made longer. (^{\wedge}\wedge;)*
 - ❖ *With increasing magnet aperture ron*
 - ❖ *For higher energy neutron*
- ❖ *Acceptable beam divergence o(^{\wedge}\wedge;)o*
 - ❖ *Increases with maximum field strength. i.e. proportional to square root of $|B|_{max}$*



Polarizing power of Quadrupole magnet with different aperture sizes

Condition : $|IB| = 2\text{T}$ on the inner surface of the magnet




The magnetic field strength at the inner surface of the magnet is set constant of $|IB|=2\text{T}$, so that the magnetic field gradient linearly decreases with the aperture size.



Characteristic of the Quadrupole Magnet Polarizer


- ❖ *High polarizing power \(^{\wedge}\circ\wedge)*
 - ❖ *High Polarization degree of $P=0.9993\pm0.0025$ was achieved.*
- ❖ *High stability o(^-\wedge)o*
 - ❖ *It is composed of the permanent magnets.*
- ❖ *High transmission d(\wedge o\wedge)b*
 - ❖ *No attenuation of polarized neutrons occurs, since they do not interact with any materials.*
- ❖ *Opening aperture is small. (T.T)*
 - ❖ *The magnet aperture can be enlarged to some extent.*
- ❖ *Magnet length must be made longer. (^^\wedge;)*
 - ❖ *With increasing magnet aperture*
 - ❖ *For higher energy neutron*



- ❖ *Acceptable beam divergence o(^_\wedge;\wedge)o*
 - ❖ *Increases with maximum field strength. i.e. proportional to square root of $|B|_{max}$*

Magnet length:
2200mm

Aperture size :
15mmΦ



Quadrupole Magnet with Extended-Halbach Magnetic Circuit


Target wavelength > 6 Å

This magnet is planed to be install into SANS-J-II of JRR-3



Characteristic of the Quadrupole Magnet Polarizer

- ❖ **High polarizing power \(^{\wedge}\circ\wedge)**
 - ❖ *High Polarization degree of $P=0.9993\pm0.0025$ was achieved.*
- ❖ **High stability o(^-\wedge)o**
 - ❖ *It is composed of the permanent magnets.*
- ❖ **High transmission d(^~o^~)b**
 - ❖ *No attenuation of polarized neutrons occurs, since they do not interact with any materials.*
- ❖ **Opening aperture is small. (T.T)**
 - ❖ *The magnet aperture can be enlarged to some extent.*
- ❖ **Magnet length must be made longer. (^^\wedge;)**
 - ❖ *For higher energy neutron*
 - ❖ *With increasing magnet aperture*



- ❖ **Acceptable beam divergence o(^_\wedge;\wedge)o**
 - ❖ *Increases with maximum field strength. i.e. proportional to square root of $|B|_{max}$*

Acceptable Beam Divergence: θ_{max}

$$\theta_{max} = \tan^{-1} \left(\frac{v_{\perp}}{v_z} \right) \approx \frac{v_{\perp}}{v_z} = \frac{\sqrt{2\mu_n |B|_{max}/m_n}}{v_z}$$

Realized with permanent magnets

$$|B|_{max} = 2 \text{ T} \Rightarrow v_{\perp} = 4.8 \text{ m/s} \Rightarrow \theta_{max} = 0.7\theta_{C,Ni}$$

$$|B|_{max} = 4 \text{ T} \Rightarrow v_{\perp} = 6.8 \text{ m/s} \Rightarrow \theta_{max} \approx \theta_{C,Ni}$$

Realized with superconducting magnets

$\theta_{C,Ni}$: Critical angle of Ni mirror guide




2) Pulsed Neutron Focusing by using Multiplet-Magnetic Lens System for TOF SANS instruments.



Intensity of Neutron Source

Application to the SANS



Conceptual design of TOF SANS Spectrometers


In order to get intensity gain by using neutron in wide wavelength range,
the total flight-path should be shortened.

Type	Spectrometer	Total Flight-Path	Sample to Detector distance, L_2
TOF	J-PARC SANS	21 m	4.5 m
	Extended-Q SANS (SNS)	18 m	4 m
Steady	SANS-J-II (JRR-3)	-	10 m
	D11(ILL)	-	40 m
	NG3(NIST)	-	13 m

The short sample to detector distance L_2 tend to deteriorate the angular resolution in low- q region for SANS spectrometer. To compensate this effect, the neutron focusing technique becomes important for the TOF SANS spectrometers.



Focusing-geometry SANS (F-SANS) to extend low-Q limit



Setup of the F-SANS instrument

Developed Magnetic Lenses



Superconducting
Sextupole Magnet



Extended Halbach
Sextupole Magnet



Sextupole Pulse
Electro-Magnet


Manufacture year	2002	2002	2002
mode	Steady	Steady	Pulse
Magnet length Z_m	2 m	2 m (0.1m/unit)	2 m
Magnetic field gradient *)	12,800 Tm ⁻² @1.8K (9,480 Tm ⁻² @4.2K)	20,100 Tm ⁻²	7,500 Tm ⁻²
Effective aperture (diameter)	ϕ 46.8 mm	ϕ 25 mm	ϕ 22 mm
Focal length Z_f **) for 6 Å	4.3 m @1.8K (3.5 m@4.2K)	3.2 m	6 m

**) Definition of the Focal length

$$Z_f = Z_m + \frac{h}{\sqrt{G\alpha m_n \lambda}} \cot\left(\frac{\sqrt{G\alpha m_n \lambda}}{h} Z_m\right)$$

Chromatic aberration!!

Parallel to Point Focusing Geometry



How to apply Magnetic Lens to Pulsed Neutrons?



1. Changing the magnetic field strength with time

- Sextupole Pulse Electro-Magnet
- Sextupole Permanent Magnet with double layer structure
- Sextupole Permanent Magnet with variable bore size

2. Multiplet Magnetic Lens System

- Change the lens function (Focusing or Defocusing) with time

Function of Sextupole Magnet




Multiplet Magnetic Lens with Spin Flippers

Definition of the Focal length


$$Z_f = Z_m + \frac{h}{\sqrt{G\alpha m_n} \lambda} \cot\left(\frac{\sqrt{G\alpha m_n} \lambda}{h} Z_m\right)$$

$$f = Z_f - Z_m \propto \sim \frac{1}{\lambda^2 Z_m}$$

Parallel to Point Focusing Geometry



Changing effective magnet length by controlling the neutron spin state.




J. Suzuki
R. Pynn

Sextupole magnet focuses positively polarized neutron and defocus negatively polarize neutron!

Experimental test

Chopper to pulse a neutron beam


Beamline : JRR-3 C3-1-2-1 NOP



Experimental setup : Double magnet system



Beamline : JRR-3 C3-1-2-1 NOP




Experimental setup : Double magnet system


Experimental Result



Experimental Result



Triple Lens System



Schematic layout of the experimental setup of the triplet lens system.

Beamlne JRR-3 C3-1-2-1 NOP




Experimental setup

Neutrons with $\lambda > \sim 7\text{\AA}$ are available at NOP beamline




Tentative Experimental Result of the Triplet Lens System




Summary

The **Multiple Lens System** is considered **very realistic, reasonable and feasible**, because

- 1) each component has been **already developed**,
- 2) the total system is relatively compact,
- 3) all components are **linearly installed into the beam line**, so that we can use pinhole and focusing geometry without making big change of experimental setup,
- 4) we can use **polarized neutrons**.





Multiple Lens System



Focusing geometry SANS Instrument with the Magnetic Lens


SANS-J-II of JRR-3




*Constructed by JAEA Soft-matter & Neutron Optics Groups
in 2005*

Demonstration at a test beamline

Beamline : C3-1-2-1 NOP of JRR-3





Focusing SANS setup
using the magnetic lens



Test of ^3He Spin Filter @ NOP beamline

T. Ino (KEK) & JAEA Neutron Optics Group (2005)



JRR-3 C3-1-2-1 NOP




Focusing SANS experiment of Ni particles by using Polarized Neutrons

JEAEA Neutron Optics Group

JRR-3 C3-1-2-1 NOP



TEM picture of Ni particles



2d Netron Intenisty $I^{\text{off}} - I^{\text{on}}$

Anisotropic scattering typical for the magnetic scattering is observed.


FAQ

Question :

How do we determine the polarization degree?

Answer:

We spatially separate positive and negative components by using the setup of Stern-Gerlach experiment. From the neutron count of each part after subtracting the BG count, We calculate the polarization degree.



Stern-Gerlach experiment


The obtained value is the product of polarization degree and analyzing power of the magnet.

Definition of the polarization degree

$$P = \frac{I_+ - I_-}{I_+ + I_-}$$



質問2b: なぜ、(普通の)スピニアナライザーを使わないのか？



$$\text{Flipping ratio } R = \frac{1 + P_i P_A (1 - D)}{1 + P_i P_A (1 - D)(1 - 2f)}$$

P_i : Polarization efficiency of the polarizer
 P_A : Polarization efficiency of the analyzer
 f : Flipping efficiency of the spin flipper
 D : Depolarization factor

How to calculate the error of P ?

P is defined as, $P = \frac{I_+ - I_-}{I_+ + I_-}$

Taking the derivatives of P , we thus find

$$\left\{ \begin{array}{l} \frac{\partial P}{\partial I_+} = \frac{1}{I_+ + I_-} - \frac{I_+ - I_-}{(I_+ + I_-)^2} = \frac{2I_-}{I_{tot}^2} \\ \frac{\partial P}{\partial I_-} = -\frac{1}{I_+ + I_-} - \frac{I_+ - I_-}{(I_+ + I_-)^2} = -\frac{2I_+}{I_{tot}^2} \end{array} \right.$$

$$\therefore I_{tot} = I_+ + I_-$$

The covariance is obviously 0 here since the measurements are independent.

$$\sigma_P^2 \cong \frac{4I_-^2}{I_{tot}^4} \sigma_{I_+}^2 + \frac{4I_+^2}{I_{tot}^4} \sigma_{I_-}^2$$

$$\sigma_P \cong \sqrt{\frac{4I_-^2}{I_{tot}^4} \sigma_{I_+}^2 + \frac{4I_+^2}{I_{tot}^4} \sigma_{I_-}^2}$$


FAQ

Question :

We are analyzing the small portion of the output beam.

What is the polarization of the whole amount of the output beam?

How can we measure the polarization of the whole amount of the output beam?




偏極度の定義

$$P = \frac{I_+ - I_-}{I_+ + I_-}$$




Polarization Analysis by using the magnetic lens



Result

Flipping ratio = $\frac{\Delta h_f}{\Delta h_{df}} = \frac{1 + P_i P_A (1 - D)}{1 + P_i P_A (1 - D)(1 - 2f)}$



Assumption : $P_A = f = 1, D = 0$


$$P_i = \frac{\Delta h_f / \Delta h_{df} - 1}{\Delta h_f / \Delta h_{df} + 1} \sim 0.999$$

Polarization Analysis by using the magnetic lens: Result 2



$$\text{Flipping ratio} = \frac{\Delta h_f}{\Delta h_{df}} = \frac{1 + P_i P_A (1 - D)}{1 + P_i P_A (1 - D)(1 - 2f)}$$

Assumption: $P_A = f = 1, D = 0$


$$P_i = \frac{\frac{\Delta h_f}{\Delta h_{df}} - 1}{\frac{\Delta h_f}{\Delta h_{df}} + 1}$$



A previous study of the Quadrupole magnet




magnetic field strength : $|BI| = G_q r$, $G = 140 \text{ T/m}$




Function of the magnets


Quadrupole



Prism



Sextupole



Lens

